Register your instrument! www.eppendorf.com/myeppendorf

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Copyright© 2013 Eppendorf AG, Hamburg. No part of this publication may be reproduced without the prior permission of the copyright owner.

Trademarks

Eppendorf®, the eppendorf logo, epMotion®, epT.I.P.S.® and are registered trademarks of Eppendorf AG.

LightCycler® and MagNA Pure® are a registered trademark of Roche Diagnostics.

Excel® and Windows® are registered trademarks of Microsoft Corporation in the United States and other countries.

Registered trademarks are not marked in all cases with ® or ™ in this manual.

The software of the device (firmware) contains open source software. License information is available on request from Eppendorf AG.

Only for epMotion with magnetic seperation: NOTICE TO PURCHASER; LIMITED LICENSE FOR RESEARCH USE ONLY.

This product and its use may be covered by one or more patents owned by Gen-Probe Incorporated. The purchase price for this product includes only limited, nontransferable rights under certain claims of certain patents owned by Gen-Probe Incorporated to use this product for research purposes only. No other rights are conveyed. Purchaser is not granted any rights under patents of Gen-Probe Incorporated to use this product for any commercial use. Further information regarding purchasing a license under patents of Gen-Probe Incorporated to use this product for any other purposes, including, without limitation, for commercial use, may be obtained by contacting Gen-Probe Incorporated, Attn: Business Development Department, 10210 Genetic Center Drive, San Diego, California 92121-4362, U.S.A./Paragraph

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Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
message without code	Make sure there is nothing in the gripper: If any labware is held by the gripper arms remove this item manually; then press Ok.	Method run was aborted; gripper is still fixed in the tool holder; before tool holder is moving further (e.g., to its home position) the message appears.	If the gripper still holds any item: 1. Remove the item manually by opening the two gripper arms manually and removing the labware manually. 2. Press button <i>Ok</i> . The gripper expects that the labware was removed and moves now to the home position.
(without code)	Because all tool holder positions are occupied the tool has to be removed manually: Press Ok, open hood and be prepared to grasp the tool manually as the tool holder opens and drops the tool.	Method run was aborted; dispensing tool is still fixed in the tool holder.	➤ To remove the dispensing tool follow the instructions.
0x0600, 0x060E	Tool did not find home	 Home position for the tool is not found. No tool inserted. Tool damaged.	 Insert tool. Check tool. Reboot and try again. If error occurs again: Call local Eppendorf Service.
0x0C01	Volume too large for this tool	Volume to be dispensed is too large for the selected tool. Possible causes: • Errors in tool files. • Errors in liquid type files.	Call local Eppendorf Application Support.
0x0C02	Volume too small for this tool	Volume to be dispensed is too small for the selected tool. Possible causes: • Errors in tool files. • Errors in liquid type files.	➤ Call local Eppendorf Application.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x0C08	Tool dimension unknown	Tool dimension values unknown. Labware outdated or corrupt.	Make sure all labware is of the latest version.
0x120A	Program aborted by user	User pressed the <i>Abort</i> button during program run.	 Error message was an information for the user that he had aborted.
0x1221	The hood was opened while the program was stopped	See error message.	➤ Close hood.
0x1222	Transfer allowance was prematurely deactivated during program initialisation	See error message.	► Start program again.
0x128B	Carrier: steps lost in x	Carrier collision. Carrier mechanics out of order.	 Ensure that the worktable equipment in epBlue is identical to the epMotion, then restart the application. If the error persists, call Eppendorf Service.
0x128C	Carrier: final position in y not found	 Problems in carrier movement in y-axis (sluggish movement or no movement at all). Light barrier for carrier in y-axis defective. 	► Call Eppendorf Service.
0x128E	Carrier: steps lost in y	Carrier collision Carrier mechanics out of order.	 Ensure that the worktable equipment in epBlue is identical to the epMotion, then restart the application. If the error persists, call Eppendorf Service.
0x1294	Carrier: steps lost in z	 Carrier collision. Carrier mechanics out of order. 	 Ensure that the worktable equipment in epBlue is identical to the epMotion, then restart the application. If the error persists, call Eppendorf Service.
0x1295	Carrier: steps lost in z before picking up tip	 Tip was still on pipet tool when tool started to pick up a new tip. Tip rack not placed correctly on the worktable. Mechanical problems of carrier. 	 Remove tips from tools. Place tip rack correctly and plane on the worktable. In other cases: Call local Eppendorf Service.
0x1296	Maximum number of tool cylces exceeded	See error message.	► Have the tool serviced.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1297	Danger of collision	 When running the programmed application, the tool carrier system will touch racks or other labware on the worktable; e.g., during pipetting the Optical Sensor may touch a long tube on the adjacent position; possible reasons: A low plate (microplate) is located next to a high tube rack. The 50 μL or 300 μL tip is programmed to move almost to the bottom of a very long tube with another long tube in the adjacent position. 	 Program the labware on the worktable in a way that high and low labware are not adjacent. Program the labware in a way that the 30 μL or 300 μL tip does not have to move deeply into a long vessel. If possible: use higher volumes in the long vessels. If possible: use longer tips for the long vessels.
0x1298	Tool not calibrated	The tool is not calibrated.	Have the tool serviced.
0x1299	Invalid number of samples	Value for Number of Samples not permissible.	Insert an admissible value for Number of Samples.
0x129A	Tip too small	Reagent Transfer: Used tip is too small.	Use a larger tip.
0x129B	Source vessel too small	Reagent Transfer: Used source vessel is too small.	Use a larger vessel.
0x12B4	Temperature unit in location C1 not active	In command wait the parameter wait for temperature was edited for location C1 although no temperature command for C1 had been edited before.	Insert command temperature before command wait or delete the parameter wait for temperature for this location.
0x12B5	Temperature unit in location C2 not active	In command wait the parameter wait for temperature was edited for location C2 although no temperature command for C2 had been edited before.	Insert command temperature before command wait or delete the parameter wait for temperature for this location.
0x12B6	Temperature unit in location C3 not active	In command wait the parameter wait for temperature was edited for location C3 although no temperature command for C3 had been edited before.	Insert command temperature before command wait or delete the parameter wait for temperature for this location.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x12B7	Temperature unit in location C1 too hot for labware (temperature command in this method)	The temperature of location C1 (according to a <i>temperature</i> command in this application) is too high for the labware placed on this location.	➤ Select a labware item that is more resistant to high temperatures (e.g. made from polypropylene or polycarbonate).
0x12B8	Temperature unit in location C2 too hot for labware (temperature command in this method)	The temperature of location C2 (according to a temperature command in this application) is too high for the labware placed on this location.	➤ Select a labware item that is more resistant to high temperatures (e.g. made from polypropylene or polycarbonate).
0x12B9	Temperature unit in location C3 too hot for labware (temperature command in this method)	The temperature of location C3 (according to a <i>temperature</i> command in this application) is too high for the labware placed on this location.	more resistant to high
0x12BA	Temperature unit in location C1 too hot for labware (temperature command in previous method)	The temperature of location C1 (according to a temperature command in a previous application with parameter keep temperature after method run set active) is too high for the labware placed on this location.	 Select a labware item that is more resistant to high temperatures (e.g. made from polypropylene or polycarbonate). Change previous application not to keep temperature after application run.
0x12BB	Temperature unit in location C2 too hot for labware (temperature command in previous method)	The temperature of location C2 (according to a temperature command in a previous application with parameter keep temperature after method run set active) is too high for the labware placed on this location.	 Select a labware item that is more resistant to high temperatures (e.g. made from polypropylene or polycarbonate). Change previous application not to keep temperature after application run.
0x12BC	Temperature unit in location C3 too hot for labware (temperature command in previous method)	The temperature of location C3 (according to a temperature command in a previous application with parameter keep temperature after method run set active) is too high for the labware placed on this location.	 Select a labware item that is more resistant to high temperatures (e.g. made from polypropylene or polycarbonate). Change previous application not to keep temperature after application run.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x12D0	Parameter conflict: Elution volume too large for this tool	Sample transfer with option elution from filter: Volume to be aspirated is too large for the tip used.	Select a tip large enough for picking up the liquid as well as the additional volume of air to be aspirated when using this option.
0x12D1	Parameter conflict: Elution volume too large for destination tube or well	Option <i>elution from filter</i> : volume is too large for the vessel used.	➤ Select a tool large enough when using this option.
0x12E1	Parameter conflict: Prewetting not possible when aspirate from bottom is selected	See error message.	► Change application.
0x12E2	Parameter conflict: Prewetting not possible when dispense from top is selected	See error message.	► Change application.
0x12E3	Parameter conflict: Prewetting not possible when elution from filter is selected	A liquid type using a prewetting step (e.g., ethanol 98%) cannot be used in combination with the parameter elution from filter in a sample transfer command.	► Change application.
0x12E4	Transport of labware: danger of collision	Danger of collision when tranporting a stack of several racks. The upper rack would collide with the carrier.	Do not pile more than 5 racks in a location.
0x12E6	Level too high	The liquid level would be higher than the vessel after dispensing.	Adjust the liquid to be dispensed to the vessel.
0x12E7	Opening the hood is not allowed when putting down tool. Please switch off power, then switch on again to restart method.	See error message.	➤ See error message.
0x12F3	Thermomixer is too hot for labware (temperature command in this method)	The Temperature unit of the thermomixer is to hot for the selected Labware	Choose a lower temperature
0x12F4	Thermomixer is too hot for labware (temperature command in previous method)	The Temperature unit of the thermomixer is to hot for the selected Labware	Choose a lower temperature
0x12F7	Waiting for thermomixer	The procedure is waiting for the thermomixer.	Wait until thermomixer function has ended

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x12F8	The selected mix speed ist not possible with this labware	The mix speed is not allowed for the selcted labware	Select another labware oder mix speed
0x1500	too big vessel index in location:	A tube is to be accessed for which the index is greater than the number of tubes on the plate/rack/holder.	► Error during creation of the application.
0x1504	<pre><rack name=""> is not accessible for tools in location</rack></pre>	Rack is a lower part of a labware stack; therefore, the tool has no access.	► Change application so that the rack is accessible.
0x1509	Liquid volume too large for vessel in location:	Total volume supplied in a source vessel is larger than needed or larger than vessel.	 Provide less volume in the vessel. Change application. When verifying the total volume needed for the source or destination, take into account additional aspirated volume in case of multidispense mode.
0x150A	Liquid volume too low for vessel in location:	Total volume supplied by the user in a source vessel is smaller than needed for a sample transfer, reagent transfer or mix command (total volume = volume to be aspirated + remaining volume for this vessel + (in case of multidispense mode:) additional aspirated volume.	 Calculate the total volume for the source or destination vessel and select a suitable vessel. Regarding additional aspirated volume in case of multidispense mode, refer to manual. Consider that the software may calculate higher remaining volumes in some cases to avoid crashes. Set Liquid detection to off for racks that are on park positions at the beginning of the procedure.
0x150B	Optical sensor: Liquid volume too low in location:	See above (error 0x150A).	➤ See above (error 0x150A).

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x150D	Optical sensor: Plate could not be found in location: Optical sensor: Rack could not be found in location:	The rack programmed for this location could not be found by the optical sensor; possible causes: Rack not placed onto location (wrong rack code or wrong rack height). Rack in wrong orientation. Problems related to the optical sensor function.	 Place the rack onto the locations as edited in the corresponding application; or: Make sure that the rack is placed plane on the worktable surface; or: Rotate rack 180° (front to back) and place it back onto the worktable location; or: Call local Eppendorf Service.
0x150E	Optical sensor: Tips could not be found in location	 The tip rack programmed for this location could not be found by the Optical Sensor; possible causes: Tip rack not placed onto location. Problems related to the Optical Sensor function. 	 Place the tip rack onto the locations as edited in the corresponding application; or: Call local Eppendorf Service.
0x1510	Optical sensor: Nothing could be found in location:	See error message.	 Place the labware programmed for this location on the worktable. If error occurs again: Call local Eppendorf Service.
0x1512	Tip type is not placed on the worktable	Tips that are needed according to the application are not available on the worktable.	 Place the tip tray programmed for this location on the worktable. If error occurs again: Call local Eppendorf Service.
0x1513	Position is out of range	The position to be addressed by the tool carrier is outside of its available range. Possible cause: Rack in park position is programmed to be addressed by the dispensing tool.	 Change application. If error occurs again: Call local Eppendorf Service.
0x1514	Optical sensor: Rack in wrong orientation in location	The tub holder has been placed onto the worktable in the wrong direction.	▶ Rotate tub holder 180° and place it back onto the worktable; restart the application.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1515	Tool cannot be used for rack in location	Distance between tip cones of the liquid handling tool does not match the distance between vessels (e.g., 24 tubes - rack does not fit the 8-channel tool).	► Change application.
0x1516	No vessel in location:	Vessels that are needed according to the application are not available on the worktable (vessel/rack combination).	 Place the vessel/rack combination programmed for this location on the worktable. If error occurs again: Call local Eppendorf Service.
0x1518	Rack is not tangible in location:	The item can not be picked up by the gripper.	► Change the application.
0x1519	Tip is too thick for vessel in location:	Diameter of the destination vessel is too small for the tip when dispensing the liquid.	 Select other tips or vessels in the application. Select dispense from top in the options of the liquid handling command.
0x151A	Optical sensor: There is a cap on vessel in location:	The Optical Sensor has detected a cap on a vessel when trying to detect a liquid level.	Remove the cap from the vessel and start the run again.
0x151B	Optical sensor: There is a wrong vessel in location:	Relates to vessels that are equipped with a readable code (e.g. Eppendorf tubs): The rack programmed for this location could not be found by the Optical Sensor; possible causes: • Wrong vessel. • Problems related to the Optical Sensor function.	 Place the vessel onto the location as edited in the corresponding application; or: Call local Eppendorf Application Support.
0x151C	Optical sensor: Vessel too high for level detection in location:	Level detection for very high vessels is not possible.	 Switch off the level detection for this vessel. Use level detection only for vessel/rack equipment with a total height below 103 mm.
0x151E	Detected volume is out of detection range	Normally a system/hardware error (malfunction of the Optical Sensor); but may also be caused by filling a vessel up to the total vessel height.	 Do not fill vessels above the specified maximum filling volume. In other cases: Call local Eppendorf Service.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x151F	Labware stack too high in location: Maximum pieces which may be piled:	A maximum of 5 racks can be stacked in a location. Placing more than 5 racks in a location.	Do not stack more than 5 racks in a location.
0x1522	Transport of the selected labware to this position or stack is not allowed Position xxx	Different reasons to forbid this transport; main reason: transport of the selected labware to this position or stack might result in instability of the stack, or the resulting stack might become too high.	► Change application.
0x1523	Labware stack too high in location: xxx Maximum height: xxx mm	Total height of labware stacks (like plates on adapters) is too high.	► Edit application to avoid too high labware stacks.
0x1526	No free position for deposition of tool available	When trying to deposit the dispensing tool after use the tool holder did not find a free position for the tool.	► Clear at least one position on the worktable to accept a dispensing tool.
0x1529	Dispensing tool could not be found. Dispensing tool required for this method:	See error message.	 Provide the dispensing tool that is needed on the worktable at its proper position.
0x152A	Gripper could not be found	See error message.	 Provide the gripper on the worktable at its proper position.
0x152D	Tip too short Select other tips or vessels in the method.	Tip does not reach the liquid level at the beginning or during the course of the liquid handling command.	Select other tips or vessels in the application.
0x152F	Labware stack too high in location: xxx Maximum pieces which may be piled: 5	A maximum of 5 racks can be piled in a location. After a <i>transport</i> command the labware stack would be too high.	► Change application.
0x1530	Dosing in pile not allowed in location: xxx	See error message.	► Change application.
0x1531	Transport of this pile not allowed in location: xxx	Do not transport a stack.	► Change application.
0x1581	Optical sensor: Liquid level could not be detected in location:	Error in level detection.	▶ Repeat measurement.
0x1602	Program too long	Maximum size of steps is reached.	 Delete commands that are not necessary before insert new commands.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1700	Liquid volume too low for vessel in location:	Total volume supplied by the user in a source vessel is smaller than needed for a sample transfer, reagent transfer or mix command (total volume = volume to be aspirated + remaining volume for this vessel + (in case of multidispense mode:) additional aspirated volume.	 Calculate the total volume for the source or destination vessel and select a suitable vessel. Regarding additional aspirated volume in case of multidispense mode, refer to manual. Consider that the software may calculate higher remaining volumes in some cases to avoid crashes. Set Liquid detection to off for racks that are on park positions at the beginning of the procedure.
0x1701	Liquid volume too large for vessel in location:	Total volume supplied in a source vessel is larger than needed or larger than vessel.	 Provide less volume in the vessel. Change application. When verifying the total volume needed for the source or destination, take into account additional aspirated volume in case of multidispense mode. epMotion 5070 only: Set Liquid detection to off for racks that are on park positions at the beginning of the procedure.
0x1904	The following labware has been deleted:	Edit mode: The worktable was changed after a application had been programmed; thus, the labware defined in a command is no longer available.	 Change the source or destination in the parameter of the respective command in accordance to match the worktable. In this case the pattern also has to be re-edited; The labware has to be re-programmed in the worktable.
0x1908	The method was written with a newer program structure. You must update your software, if you want to edit this method	See Error message.	Update your software, or:Call local Eppendorf Service.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x190D	The following labware is not selected in the Labware File Window:	Edit mode/worktable: The chosen labware is not available in the labware collection that had been selected for your lab. Possible cause for this error message: The labware has been de-selected in the Labware File Window.	▶ Select the respective labware in the Labware File Window. You need to have the concering User rights. If you don't have the needed User rights ask your administrator.
0×190E	The following tool is not selected in the Labware File Window:	Edit mode/procedure: The chosen tool is not available in the labware collection that had been selected for your lab. Possible cause for this error message: The labware has been de-selected in the Labware File Window.	➤ Select the respective labware in the Labware File Window. You need to have the concering User rights. If you don't have the needed User rights ask your administrator.
0x1910	The position of the following labware is not available on this worktable The method was written for another workstation configuration	Edit mode/worktable: The chosen labware is not available in the labware collection that had been selected for your lab in the software node <i>epMotion</i> . Possible cause for this error message: The labware has been de-selected in the <i>epMotion</i> node.	➤ Set the respective labware in the node <i>epMotion</i> .
0x1910	The method was written for another workstation configuration The position of the following labware is not available on this worktable	The position of required labware is not available on this device e.g. you have a 5075 LH and the application was written on a 5075 VAC.	 Load the concerned application on an compatible device, or Modify the application until it fits for the available device.
	The position of the following labware is not allowed anymore yyy in location: xxx	Edit mode/worktable: The chosen labware may not be placed on the selected position anymore. Possible cause for this error message: The application has been written with a former version of the software.	Place the respective labware on another position.
0x1911	The following labware has been changed, so that the pattern does not fit anymore: xxx	Because you can change the order or contents <i>Tubs</i> + <i>Modules</i> (<i>equip</i>) + <i>Holders</i> combination, it can happen that the recent pattern of a command does not fit the new positions of the tubes.	▶ Either change the order or contents of the <i>Tubs</i> + <i>Modules(equip)</i> + <i>Holders</i> combination back to the original. Or change the pattern in the command.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1911	The following labware has been changed, so that the pattern does not fit anymore: xxx	As it is possible to change the order or contents of a "Equiped Holder"-combination, it can happen that the recent pattern of a command does not fit the new positions of the tubes.	▶ Either change the order or contents of the "Equiped Holder"-combination back to the original. Or change the pattern in the command.
0x1A02	The name is already used for another labware	Edit mode / labware: The same name has been defined for a different rack or another labware item.	► Enter a different name.
0x1A03	This position is not available for the selected labware	Edit mode/worktable: Certain worktable positions are not allowed for certain labware (e.g., tips can only be placed in the rear of the worktable).	▶ Place the selected labware in another location.
0x1A04	The selected labware may not be stacked on top of labware already placed	Edit mode/worktable: Building of labware stacks on the worktable is restricted to certain labware combinations (e.g., thermorack above thermorack does not make sense).	➤ See "Cause".
0x1A05	Maximum number of labware stacked: 5	Edit mode/worktable: Certain labware build a pile which contains too many when placed on certain labware.	Place the selected labware in another location.
0x1A06	Labware stack too high in location: xxx Maximum height: xxx mm	Edit mode/worktable: Labware stacks on the worktable may not exceed a maximum height limit (e.g., plates on adapters is allowed; reservoir holder on adapters is not allowed because the stack would become too high).	➤ See "Cause".
0x1A10	8 channel tool cannot be used for this source rack	Edit mode/parameter in command Sample Transfer: Source rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	► Choose another rack or another tool.
0x1A10	8 channel tool cannot be used for this source rack	Edit mode / parameter in command Sample Transfer: Source rack does not fit to 8-channel-tool (e.g.: 24-well-plate or tube rack with 24 positions).	Choose another rack or another tool.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1A11	8 channel tool cannot be used for this destination rack.	Edit mode / parameter in command Sample Transfer: Destination rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	➤ Choose another rack or another tool.
0x1A12	No source or destination selected	Edit mode/parameter in command <i>Sample Transfer</i> : Source or destination rack has not been selected.	Select source or destination, respectively.
0x1A15	Invalid entry for movement blow (0 100)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>Movement Blow</i> .	► Enter a value between 0 and 100%.
0x1A15	Invalid entry for movement blow (0 100)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter movement blow.	► Enter a value between 0 and 100%.
0x1A16	Invalid entry for delay blow (0 9999)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>Delay Blow</i> .	► Enter a value between 0 and 9999 msec.
0x1A16	Invalid entry for delay blow (0 9999)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>delay blow</i> .	► Enter a value between 0 and 9999 msec.
0x1A17	Invalid entry for speed aspiration (0.2 110)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>Speed Aspiration</i> .	► Enter a value between 0.2 and 110 mm/sec.
0x1A17	Invalid entry for speed aspiration (0.2 110)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter speed aspiration.	► Enter a value between 0.2 and 110 mm/sec.
0x1A18	Invalid entry for speed dispense (0.2 110)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>Speed Dispense</i> .	► Enter a value between 0.2 and 110 mm/sec.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1A18	Invalid entry for speed dispense (0.2 110)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter speed dispense.	► Enter a value between 0.2 and 110 mm/sec.
0x1A19	Invalid entry for speed blow (0.2 110)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>Speed Blow</i> .	► Enter a value between 0.2 and 110 mm/sec.
0x1A19	Invalid entry for speed blow (0.2 110)	Edit mode/parameter in transfer command: A value beyond the allowed range has been entered for the parameter <i>speed blow</i> .	► Enter a value between 0.2 and 110 mm/sec.
0x1A1A	Invalid entry for initial stroke (0 100)	Edit mode/parameter in transfer command: a value beyond the allowed range has been entered for the parameter initial stroke.	► Enter a value between 0 and 100 %.
0x1A20	8 channel tool cannot be used for this source rack	Edit mode/parameter in command <i>Reagent Transfer</i> : Source rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	Choose another rack or another tool.
0x1A21, 0x1A30, 0x1A31, 0x1A51	8 channel tool cannot be used for this destination rack	Edit mode/parameter the affected command: Destination rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	► Choose another rack or another tool.
0x1A22, 0x1A32, 0x1A42, 0x1A52	No source or destination selected	Edit mode/parameter in the affected command: Source or destination rack has not been selected.	Select source or destination, respectively.
0x1A40, 0x1A41, 0x1A50	8 channel tool cannot be used for this source rack	Edit mode/parameter in the affected command: Source rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	Choose another rack or another tool.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1A61	8 channel tool cannot be used for this rack	Edit mode/parameter in command <i>Mix</i> : Rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	➤ Choose another rack or another tool.
0x1A62	No rack selected	Edit mode/parameter in command <i>Mix</i> : Rack has not been selected.	► Select rack, respectively.
0x1A65	Invalid entry for speed (0.2 110)	Edit mode/parameter in command <i>Mix</i> : A value beyond the allowed range has been entered for the parameter <i>speed</i> .	Enter a value between 0.2 and 110 mm/sec.
0×1A70	This position is already occupied	Edit mode/pattern: When editing the pattern you have tried to select a certain position that is already occupied.	Follow the direction of the edited pattern and move to a diffferent position.
0x1A76	8 channel tool cannot be used for this module rack	Edit mode/pattern: Rack does not fit the 8-channel tool (e.g., <i>Tubs</i> + <i>Modules</i> (equip) + <i>Holders</i> -combination with positions all less than 8 in modules).	➤ Choose another rack or another tool.
0x1A77	No module rack or tubes found	Edit mode/pattern: Rack does not have any positions (e.g., <i>Tubs + Modules</i> (equip) + Holders-combination with positions all less than 1 in Modules).	► Choose another rack.
0x1A78	Number of tubes not supported	Edit mode/pattern: One or more Modules have 3, 5, 6, 7 or more than 8 positions. This is not supported.	► Choose another rack.
0x1A80	Invalid entry for minutes (0 99)	Edit mode/parameter in command <i>Wait</i> : A value beyond the allowed range has been entered for the parameter <i>minutes</i> .	➤ Enter a value between 0 and 99 minutes.
0x1A81	Invalid entry for seconds (0 59)	Edit mode/parameter in command <i>Wait</i> : A value beyond the allowed range has been entered for the parameter <i>seconds</i> .	► Enter a value between 0 and 59 seconds.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1A90	Selecting more than one rack as Source or as Destination: All source racks (or all destination racks, resp.) must have the same well pattern	Edit mode/parameter in transfer command: If more than one rack is selected as source (or destination) all source (or destination), racks must have the same well pattern (e.g. 96-well plates).	➤ See explanation in "Cause".
0x1A90	Selecting more than one rack as Source or as Destination: All source racks (or all destination racks, resp.) must have the same well pattern	Edit mode/parameter in transfer command: The selected labwares have a different amount of wells.	► Choose another rack.
0x1AA0	Invalid entry for temperature (0 110)	Edit mode/parameter in command <i>Temperature</i> : A value beyond the allowed range has been entered for the parameter <i>Temperature</i> .	► Enter a value between 0 and 110 degrees.
0x1AE0	Invalid entry for mix speed (0 2000)	Edit mode / parameter in command "Thermomixer": a value beyond the allowed range has been entered for the parameter "speed"	Enter a value between 0 and 2000 rpm
0x1AE1	Invalid entry for minutes (0 120)	Edit mode / parameter in command "Thermomixer": a value beyond the allowed range has been entered for the parameter "minutes"	Enter a value between 0 and 120 minutes
0x1AE2	Invalid entry for seconds (0 59)	Edit mode / parameter in command "Thermomixer": a value beyond the allowed range has been entered for the parameter "seconds"	Enter a value between 0 and 59 seconds
0x1AE3	Invalid entry for temperature (4 110)	Edit mode / parameter in command "Thermomixer": a value beyond the allowed range has been entered for the parameter "temperature"	Enter a value between 4 and 110 degrees

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1C0B	Sample number too large	Run mode: The number of samples you entered will fill more than one rack (source or destination, respectively) based on the programmed pattern.	 Start the application again and enter a lower number of samples; or: Enter the edit mode and program a pattern that together with the number of samples you want to run will not extend beyond one rack.
0x1C0D	You must clear old pattern first Press new pattern	Edit mode/pattern: You tried to change a stored pattern before deleting the old pattern.	▶ Delete the old pattern by pressing the softkey <i>New</i> pat.
0x1C0D	You must clear old pattern first Press "new pattern"	Edit mode/pattern: You tried to change a stored pattern before deleting the old pattern.	▶ Delete the old pattern by pressing the button new pattern.
0x1C0E	You must go forward	Edit mode/pattern: When entering the pattern, the order of edited locations in the source (or destination, respectively) must be from left to right or from the top of the pattern downwards (i.e., move only in columns or in rows).	➤ See explanation in "Cause".
0x1C0F	You may only move horizontally or vertically	Edit mode/pattern: When entering the pattern the order of edited locations in the source (or destination, resp.) must be from the left to right or from top of the pattern downwards (i.e., move only in columns or in rows). Note: Error message may also occur when working with an 8-channel tool and editing another position than the upper ones (see error code 0x1C1F).	
0x1C10	Pattern for replicates of first sample too complex		 Enter a simpler pattern if possible. In case this is not possible: Call local Eppendorf Application Support.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1C11	Pattern too complex	Edit mode/pattern: The pattern algorithm cannot handle this pattern. Note: See note in error 0x1C0F.	 Enter a simpler pattern if possible. In case this is not possible: Call local Eppendorf Application Support.
0x1C14	Pattern must fit in rows or columns	Edit mode/pattern: The basic unit of the pattern you tried to enter extends beyond a row or a column. This cannot be handled by the pattern algorithm.	 Enter a simpler pattern if possible. In case this is not possible: Call local Eppendorf Application Support.
0x1C15	Pattern too complex	Edit mode/pattern: The pattern algorithm cannot handle this pattern. Note: See note in error 0x1C0F.	 Enter a simpler pattern if possible. In case this is not possible: Call local Eppendorf Service.
0x1C16	This position is already occupied	Edit mode / pattern: When editing the pattern you have tried to select a certain rack position that is already occupied.	➤ Following the edited pattern move to a different rack position.
0x1C17	You must start with the source	Edit mode/pattern: When editing a pattern you must start with the source.	► See explanation in "Cause".
0x1C18	Please enter a source now	Edit mode/pattern: In the destination rack, you tried to enter more replicates than you had sources.	► Enter the same number of replicates for all sources you edit when programming a pattern.
0x1C19	Please enter a destination now	Edit mode/pattern: When having selected a source in the command <i>Sample Transfer</i> you first have to enter a destination for this source before moving to the next source position.	Press softkey Destination and edit the destination position(s) for the selected source position.
0x1C1A	No more positions available (limited by Number of Samples Command)	Edit mode/pattern: Editing further positions is not possible because the limit set in the command <i>Number of Samples</i> would be exceeded.	➤ Select a pattern that fits the programmed command Number of Samples.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1C1E	Pattern for reagent transfer: source can be chosen only once	Edit mode/pattern for command <i>Reagent Transfer</i> : After having entered the source and the destinations for the reagent transfer you cannot select an additional source.	 Enter the source only once. In case this does not meet your goals for this application consider selecting command Sample Transfer instead of Reagent Transfer; or: Call local Eppendorf Application Support.
0x1C1F	Pattern with 8-channel tool: Please edit upper position of this tool	Edit mode/pattern with 8-channel tool: Only the upper positions of the 8-channel tools can be selected.	➤ See explanation in "Cause".
0x1C20	Pattern for sample transfer: only one position per sample on source	Edit mode/pattern for command <i>Sample Transfer</i> : Before selecting a second source position, you have to edit the destination for the first source position.	▶ Enter destination for the source you just selected; afterwards, you can edit the next source position.
0x1C21	In source rack further positions cannot be edited because positions in destination rack are already occupied	Edit mode/pattern: Selecting further source positions would require a second destination rack according to the pattern you edited.	▶ Edit a pattern that does not require more than one destination rack per command. To use more destination racks, create additional commands.
0x1C22	Pattern for pool one dest: destination can be chosen only once	Edit mode/pattern for command <i>PoolOneDest</i> : After having entered the sources and the destination, you cannot select an additional destination.	 Enter the destination only once. In case this does not meet your goals for this application, consider selecting command <i>Pool</i> instead of <i>PoolOneDest</i>; or: Call local Eppendorf Application Support.
0x1C23	Pattern for dilute: only one position per sample on source	Edit mode/pattern for command <i>Dilute</i> : Before selecting a second source position, you have to edit the destination for the first source position.	▶ Enter destination for the source you just selected; afterwards, you can edit the next source position.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x1C25	Pattern for pool: only one position per sample on destination	Edit mode/pattern for command <i>Pool</i> : Before selecting a second destination position, you have to edit the next source positions to be pooled into this destination.	▶ Enter sources for the next destination position; afterwards, you can edit the next destination position.
0x1C26	Pattern for <i>Reagent Transfer</i> : not enough source positions	Run mode: To provide enough reagent volume for the number of samples you entered, the selected reagent source positions must be higher.	 Start the application again and enter a lower number of samples; or: Enter the edit mode and program more reagent source positons in the pattern. Keep in mind that the selected reagent source positions may not extend beyond one rack.
0x2025	Bottom tolerance too big	Bottom tolerance too big.	▶ Use a smaller value.
0x2026	Bottom tolerance too small	Bottom tolerance too small.	Use a bigger value.
0x2101	Tool not defined.	Parameter <i>Pipet. Tool</i> was not edited in the application.	➤ See explanation in "Cause".
0x2104	Tips not edited in worktable/procedure	Tips were edited in the procedure of the application, but they were not edited in the worktable (e.g., filter tips <-> tips without filter).	► Edit the tips that you programmed in the procedure in the worktable.
0x2105	Parameter conflict: Startvolume greater than filling volume of source tube or well	The parameters <i>Volume</i> and <i>Source</i> of the source vessel do not match (<i>Volume</i> is higher than the maximum filling volume of the source vessel).	▶ Edit Volume and Source in the worktable so that Volume is covered by the maximum filling volume of the source vessel.
0x2106	Parameter conflict: Startvolume greater than filling volume of destination tube or well	The parameters <i>Volume</i> and <i>Destination</i> of the destination vessel do not match (<i>Volume</i> is higher than the maximum filling volume of the destination vessel).	▶ Edit Volume and Destination in the worktable so that Volume is covered by the maximum filling volume of the destination vessel.
0x2107	Volume not defined	Parameter <i>Volume</i> was not edited in the application.	► See explanation in "Cause".
0x2108	Parameter conflict: Volume too small for this tool	The parameters <i>Volume</i> and <i>Pipet</i> . <i>Tool</i> of the application do not match (<i>Volume</i> is smaller than the lower limit of the tool volume range).	▶ Edit Volume and Pipet. Tool so that Volume is covered by the volume range of the pipet tool.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x2109	Parameter conflict: Volume too large for this tool	The parameters <i>Volume</i> and <i>Pipet. Tool</i> of the application do not match (<i>Volume</i> is higher than the upper limit of the tool volume range).	▶ Edit Volume and Pipet. Tool so that Volume is covered by the volume range of the pipet tool.
0x210A	Parameter conflict: volume greater than filling volume of source tube or well	The parameters <i>Volume</i> and <i>Source</i> of the application do not match (<i>Volume</i> is higher than the maximum filling volume of the source vessel).	▶ Edit Volume and Source so that Volume is covered by the maximum filling volume of the source vessel.
0x210B	Parameter conflict: volume greater than filling volume of destination tube or well	The parameters <i>Volume</i> and <i>Destination</i> of the application do not match (<i>Volume</i> is higher than the maximum filling volume of the destination vessel).	▶ Edit Volume and Destination so that Volume is covered by the maximum filling volume of the destination vessel.
0x210D	Source rack not defined	Parameter <i>Source</i> was not edited in the application.	► See explanation in "Cause".
0x210E	Source rack not edited in worktable	The source rack you edited in the procedure of the application has been removed from the worktable.	▶ Edit the rack that you programmed in the procedure as <i>Source</i> in the worktable, or edit a different source rack in the application.
0x210F	Source rack not selected in the Labware File Window	The source rack you edited in the application is not selected or removed in the Labware File Window and therefore is not available for programming.	➤ Select the rack in the Labware File Window or edit a different rack in the application.
0x2110	Destination rack not defined	Parameter <i>Source</i> was not edited in the application.	► See explanation in "Cause".
0x2111	Destination rack not edited in worktable	Destination rack was edited in the procedure of the application, but it was not edited in the worktable.	▶ Edit the rack that you programmed in the procedure as <i>Destination</i> in the worktable or edit a different destination rack in the application.
0x2112	Destination rack not selected in the Labware File Window	The destination rack you edited in the application is not selected or removed in the Labware File Window and therefore is not available for programming.	➤ Select the rack in the Labware File Window or edit a different rack in the application.
0x2113	Pattern not defined	Parameter <i>Pattern</i> was not edited in the application.	► See explanation in "Cause".

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x2114	Loading error (Invalid entry in pattern)	Normally a system error; but may also be caused by editing a pattern without destination positions; or: File damaged.	 Edit a pattern with source and destination positions. In other cases: Call local Eppendorf Service.
0x211B	Liquid type not defined	Parameter <i>Liquid Type</i> was not edited in the application.	Choose a "Liquid Type" for the Liquid Handling Command.
0x211D	Mixing cycles in source not defined	Parameter <i>No. of Cycles</i> in a mix procedure was not edited for the source in the application (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	➤ See explanation in "Cause".
0x211E	Invalid entry for mixing cycles in source (1 99)	Entry for the parameter <i>No. of Cycles</i> in a mix procedure for source vessels was higher than the max. limit (1 up to 99 cycles) (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	► Enter a number between 1 and 99 for the <i>No. of Cycles</i> parameter.
0x211F	Invalid entry for mixing speed in source (1 10)	Entry for the parameter <i>Speed</i> in a mix procedure for source vessels was higher than the max. limit (1 up to 10) (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	► Enter a number between 1 and 10 for the <i>Speed</i> parameter.
0x2120	Mixing volume in source not defined	Parameter <i>Volume</i> in a mix procedure for source vessels was not edited in the application (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	➤ See explanation in "Cause".

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x2121	Parameter conflict: mixing volume in source too large for this tool	The parameters <i>Volume</i> and <i>Pipet</i> . <i>Tool</i> of a mix procedure for source vessels are not in agreement (<i>Volume</i> is higher than the upper limit of the tool's volume range) (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	▶ Edit Volume and Pipet. Tool so that Volume is within pipet tool's volume range.
0x2122	Parameter conflict: mixing volume in source too small for this tool	The parameters <i>Volume</i> and <i>Pipet. Tool</i> of a mix procedure for source vessels do not match (<i>Volume</i> is less than the lower limit of the tool's volume range) (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	▶ Edit <i>Volume</i> and <i>Pipet. Tool</i> so that <i>Volume</i> is within the pipet tool's volume range.
0x2123	Parameter conflict: mixing volume in source greater than filling volume of source tube or well	The parameters <i>Volume</i> and <i>Source</i> of a mix procedure in the application do not match (<i>Volume</i> is higher than the maximum filling volume of the source vessel) (mix procedure as defined in a command <i>Mix</i> or as part of a liquid transfer command via parameter <i>Options</i>).	▶ Edit Volume and Source so that Volume is within the allowable filling volume of the source vessel.
0x2124	Mixing cycles in destination not defined	Parameter No. of Cycles in a mix procedure for destination vessels was not edited in the application (mix procedure as part of a liquid transfer command via parameter Options)	➤ See explanation in "Cause".
0x2125	Invalid entry for mixing cycles in destination (1 99)	Entry for the parameter No. of Cycles in a mix procedure for destination vessels was higher than the max. limit (1 up to 99 cycles) (mix procedure as part of a liquid transfer command via parameter Options).	► Enter a number between 1 and 99 for the <i>No. of Cycles</i> parameter.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x2126	Invalid entry for mixing speed in destination (1 10)	Entry for the parameter <i>Speed</i> in a mix procedure for destination vessels was higher than the max. limit (1 up to 10) (mix procedure as part of a liquid transfer command via parameter <i>Options</i>).	► Enter a number between 1 and 10 for the <i>Speed</i> parameter.
0x2127	Mixing volume in destination not defined	Parameter <i>Volume</i> in a mix procedure for destination vessels was not edited in the application (mix procedure as part of a liquid transfer command via parameter <i>Options</i>).	➤ See explanation in "Cause".
0x2128	Parameter conflict: mixing volume in destination too large for this tool	The parameters <i>Volume</i> and <i>Pipet</i> . <i>Tool</i> in a mix procedure for destination vessels do not match (<i>Volume</i> is higher than the upper limit of the tool volume range) (mix procedure as part of a liquid transfer command via parameter <i>Options</i>).	▶ Edit Volume and Pipet. Tool so that Volume is within the pipet tool's volume range.
0x2129	Parameter conflict: mixing volume in destination too small for this tool	The parameters <i>Volume</i> and <i>Pipet</i> . <i>Tool</i> in a mix procedure for destination vessels do not match (<i>Volume</i> is lower than the minimum allowed volume) (mix procedure as part of a liquid transfer command via parameter <i>Options</i>).	▶ Edit <i>Volume</i> and <i>Pipet. Tool</i> so that <i>Volume</i> is within the volume range of the pipet tool.
0x212A	Parameter conflict: mixing volume in destination greater than filling volume of destination tube or well	The parameters <i>Volume</i> and <i>Destination</i> in a mix procedure for destination vessels do not match (<i>Volume</i> is higher than the maximum filling volume of the destination vessel) (mix procedure as part of a liquid transfer command via parameter <i>Options</i>).	▶ Edit Volume and Destination so that Volume is within the maximum filling volume of the destination vessel.

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x212C	Parameter conflict: mix after dispense not allowed in multidispense mode	When the parameter transfer type is set to multidispense, the parameter mix after dispense cannot be edited for this command.	 Change parameter transfer type to pipette; or: Omit the mixing step; in this case you could also edit another mixing step as a new command (Mix), which would be performed after the previous command of the procedure has ended.
0x212D	Parameter conflict: 8 channel tool cannot be used for this source rack	Edit mode / parameter in liquid handling command: Source rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	➤ Choose another rack or another tool.
0x212E	Parameter conflict: 8 channel tool cannot be used for this destination rack	Edit mode/parameter in liquid handling command: Destination rack does not fit the 8-channel tool (e.g., 24-well plate or tube rack with 24 positions).	➤ Choose another rack or another tool.
0x212F	Sample number too large	The number of samples you have entered will fill more than one rack (source or destination, respectively) based on the programmed pattern.	▶ Program a pattern that together with the number of samples you want to run will not extend beyond one rack, or choose a lower number of samples in the corresponding Number of Samples command.
0x2130	Parameter conflict: tip cannot be used for this source rack	The parameter Source rack in the liquid handling command does not match the selected tool (e.g., 384-well plate and TS_1000 or TM1000_8).	Choose another rack or another tool.
0x2131	Parameter conflict: tip cannot be used for this destination rack	The parameter <i>Destination rack</i> in the liquid handling command does not match the selected tool (e.g., 384-well plate and TS_1000 or TM1000_8).	➤ Choose another rack or another tool.
0x2132	Invalid number of samples (1 384)	Number of samples you have entered is too high.	► Enter a maximum number of samples up to 384.
0x2133	Rack or plate not defined	Edit mode/parameter in command <i>Transport</i> : Labware to be transported was not edited.	➤ See explanation in "Cause".

Code	Symptom/Meldung	Mögliche Ursache	Abhilfe
0x2134	Rack not selected in the Labware File Window	The rack you edited in a command <i>Transport</i> (parameter <i>labware</i>) is not selected or removed in the Labware File Window and therefore is not available for programming.	➤ Select the rack in the Labware File Window or edit a different rack in the application.
0x2135	Location for destination not defined	Edit mode/parameter in command <i>Transport</i> : The destination location for the labware to be transported was not edited.	➤ See explanation in "Cause".
0x2136	Invalid entry for seconds (1 59)	Edit mode/parameter in command <i>Wait</i> : A value beyond the allowed range has been entered for the <i>seconds</i> parameter.	► Enter a value between 0 and 59 seconds.
0x2137	Invalid entry for minutes (1 99)	Edit mode/parameter in command <i>Wait</i> : A value above the allowable maximum minutes has been entered for the <i>minutes</i> parameter.	• Enter a value between 0 and 99 minutes.
0x2138	Method without active commands	application contains only passive commands (like wait, comment, etc.).	Insert at least one active command.
0x2139	Parameter conflict: mix before aspirating not allowed in multiaspirate mode	Pool/Pool One Destination: Mix before aspirating not allowed in multiaspirate mode.	▶ Do not mix.
0x2170	Parameter conflict: Parameter elution from filter is only possible when filter plates have been selected as source	To edit the option <i>elution from filter</i> a filter plate must have been edited as source.	➤ See explanations in "Cause".
0x2171	Parameter conflict: Multidispense mode is not allowed when selected elution from filter	See explanation in the error message.	➤ See explanation in the error message.
0x2172	Parameter conflict: Transfer volume must be set to zero when Parameter elution from filter has been selected	Using the option elution from filter, the complete volume contained in the filter plate wells is always aspirated; therefore, editing a volume to be transferred is not possible.	➤ Set the volume to zero because the entry will not have an effect in the application run.

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